TECHNICAL DATA

Glatt Type GPCG 1.1
"Glatt Powder Coater Granulator"
Fluidbed Dryer, Coater & Spray Granulator

This mobile Laboratory unit is designed to handle batches of between 0.5 and 5 litres (1 kg to 3 kg) and is suitable for drying, moistening, coating and spray-granulation.

Product to be dried is fed into the stainless-steel product container and placed onto the base cabinet. The container is inserted into the drying column and the whole assembly is hermetically sealed by a pneumatic ram, which pushes the individual parts up against their respective seals.

Air can now be drawn through the assembly, drying, fluidising and blending the product. The airflow is produced by a 2.2 kW fan – which is passed through a 4kW heating register where it is heated to the required temperature.

The heated air then flows through the perforated base of the product container and through the product. From time to time, the airflow is interrupted and a pneumatic shaking device which shakes the fabric filter to ensure that any product which may be lodged in the filter is returned to the product container.

This Plant is designed for local exhaust ventilation. A sight glass and light source allows observation of the process.

When the drying process is finished, the product container is removed from the filter housing by lowering and swinging it away where it can be tipped for emptying and further processing.

The machine is designed suitable for Spray Granulation and Coating (HS Wurster) and currently includes the following:

- 1 x 5.0 L Container for Drying and Top Spray Granulation including spray lance and nozzle.
- 1 x 5.0 L Container for Drying and Bottom Spray 'HS Wurster' Coating including spray lance and nozzle.
- 1 x 2.5 L Container for Drying and Bottom Spray 'HS Wurster' Coating including spray lance and nozzle.
- 1 x Watson-Marlow 520U Peristaltic Pump
- 1 x Mettler-Toledo XS1003S Precision Balance.
- 1 x Carrying case with various accessories.
- 1 x Rolling trolley for accessories

Operator controls and electronic indicators are mounted in a Glatt Ecoview II® microprocessor control housed in a Siemens SIMATIC OP37 panel to the right-hand side.

All functions are programmed, adjusted, controlled and monitored via the control. The operating language can be switched between English and German.

All drives, fans, heating elements, filters, pneumatic and electric installations are housed within the machine base.

Connection to a supply of compressed air is required.

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